INVENTUM

maakt 't moment

bloeddrukmeter • bloodpressure monitor

BDA432

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English

1	safety instructions	page 32
2	appliance description	page 37
3	things to know about blood pressure	page 39
4	before first use	page 40
5	measuring blood pressure	page 44
6	evaluating results	page 45
7	saved records	page 46
8	delete all saved records	page 48
9	tips for measurement	page 49
10	cleaning & maintenance	page 49
11	frequently asked questions	page 50
12	troubleshooting	page 50
13	specifications	page 52
14	complied European standards list	page 53
15	EMC guidance	page 55
	<i>Guarantee & after sales Service</i>	page 63

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1 safety instructions

The signs below might be in the user manual, on the labelling or on the appliance. The description of the symbols can be found in the table below.

	Refer to instruction manual/booklet To signify that the instruction manual/ booklet must be read.		Type BF applied part
CE 0123	CE marking indicates that a product has been assessed by the manufacturer and deemed to meet EU safety, health and environmental protection requirements. It is required for products manufactured anywhere in the world that are then marketed in the EU.		Caution Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences
	Symbol for "MANUFACTURER"		Symbol for "Authorised Representative in the European Community"
	The Green Dot is the license symbol of a European network of industry-funded systems for recycling the packaging materials of consumer goods.		The symbol indicates that the product should not be discarded as unsorted waste but must be sent to separate collection facilities for recovery and recycling.
	Symbol for "MANUFACTURE DATE"		Symbol for "DIRECT CURRENT"

CAUTION

- This device is intended for adult use only.
- This appliance is intended for indoor use only.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure.
- It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- This device is intended for the indication of the blood pressure and heartbeat only. A test for clinical results can be done by a physician.
- Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Do not begin or end medical treatment without asking a physician for treatment advice.
- If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your physician.
- When the device was used to measure patients who have common arrhythmias, such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur with a deviation. Please consult your physician about the result.
- When the cuff pressure reaches 40 kPa (300 mmHg), the unit will automatically deflate. If the cuff not deflate when pressures reach 40 kPa (300 mmHg), detach the cuff from the arm and press the START/STOP button to stop inflation.
- The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anaesthetic mixture with air of oxygen or nitrous oxide.
- The operator shall not touch the output of batteries and the patient simultaneously.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- This device is contraindicated for any female who may be suspected of being, or is, pregnant. Besides providing inaccurate readings, the effects of this device on the foetus are unknown.
- This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become numb, swollen and even purple due to a lack of blood.
- Please use the device in an environment which is indicated in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.

- During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential allergic reaction or contact injury.
- Please use the ACCESSORIES and detachable parts specified/authorised by the MANUFACTURER. Otherwise, it may cause damage to the unit or danger to the user/patients.
- The device doesn't need to be calibrated within two years of reliable service.
- Please dispose of the ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- To avoid measurement errors, carefully read this manual before using the product.
- Please use the soft cloth to clean the whole unit. Don't use any abrasive or volatile cleaners.
- The device is not intended for patient transport outside a healthcare facility.
- Don't kink the connection tube during use, otherwise, the cuff pressure may continuously increase which can prevent blood flow and result in harmful injury to the PATIENT.
- When using this device, please pay attention to the following situation which may interrupt blood flow and influence blood circulation of the patient, thus cause harmful injury to the patient: connection tubing kinking too frequent and consecutive multiple measurements; the application of the cuff and its pressurization on any arm where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present; inflating the cuff on the side of a mastectomy.
- Warning: Do not apply the cuff over a wound; otherwise it can cause further injury.
- Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, because this could cause temporary loss of function of those simultaneously-used monitoring ME equipment.
- On the rare occasion of a fault causing the cuff to remain fully inflated during measurement, open the cuff immediately. Prolonged high pressure (cuff pressure > 300mmHg or constant pressure > 15mmHg for more than 3 minutes) applied to the arm may lead to an ecchymosis.
- Please check that operation of the device does not result in prolonged impairment of patient blood circulation.
- When measurement, please avoid compression or restriction of the connection tubing.
- The device cannot be used with HF surgical equipment at the same time.
- The ACCOMPANYING DOCUMENT shall disclose that the SPHYGMOMANOMETER was clinically investigated according to the requirements of ISO 81060-2:2013.
- To verify the calibration of the AUTOMATED SPHYGMOMANOMETER, please contact the manufacturer.
- This device is contraindicated for any female who may be suspected of, or is pregnant.

Besides providing inaccurate readings, the effects of this device on the fetus are unknown.


- Too frequent and consecutive measurements could cause disturbances in blood circulation and injuries.
- When not in use, store the device in a dry room and protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on the storage case.
- This device may be used only for the purpose described in this booklet. The manufacturer cannot be held liable for damage caused by incorrect application.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in this booklet.
- The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Warning: No servicing/maintenance while the ME equipment is in use.
- The patient is an intended operator.
- The patient can measure data and change batteries under normal circumstances and maintain the device and its accessories according to the user manual.
- The blood pressure monitor and the cuff are suitable for use within the patient environment. If you are allergic to polyester, nylon or plastic, please don't use this device.
- If you experience discomfort during a measurement, such as pain in the arm or other complaints, press the START/STOP button to release the air immediately from the cuff. Loosen the cuff and remove it from your arm.
- Before use, make sure the device functions safely and is in proper working condition. Check the device, do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Keep the unit out of reach of infants, young children or pets to avoid inhalation or swallowing of small parts. It is dangerous or even fatal.
- Be careful to strangulation due to cables and hoses, particularly due to excessive length.
- At least 30 min required for ME equipment to warm from the minimum storage temperature between uses until it is ready for intended use. At least 30 min required for ME equipment to cool from the maximum storage temperature between uses until it is ready for intended use.
- This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;

- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d away from the equipment. The distance d is calculated by the MANUFACTURER from the 80 MHz to 5.8 GHz column of Table 4 and Table 9 of IEC 60601-1-2:2014, as appropriate.
- Please use ACCESSORIES and detachable parts specified/ authorized by MANUFACTURE. Otherwise, it may cause damage to the unit or danger to the user/patients.

2 appliance description

Component list of pressure measuring system

- 1 Cuff
- 2 Air pipe
- 3 PCBA
- 4 Pump
- 5 Valve


The packaging has the following content:

- Blood pressure monitor
- Cuff (Type BF applied part) 22-32 cm
- 4 AAA-batteries
- User manual
- Storage pouch

A larger cuff (22-42 cm) can be ordered separately. Check the website: www.inventum.eu for more information.

The following symbols can be found on the display:


Symbol	Description	Explanation
SYS	Systolic blood pressure	High pressure result
DIA	Diastolic blood pressure	Low pressure result
	Pulse	Pulse/minute: heartbeats/minute
	Deflating	Cuff is deflating
	Memory	Shows the saved measurements in the memory
	Low battery	Batteries are low and need to be replaced
kPa	kPa	Measurement unit of the blood pressure(1 kPa=7.5 mmHg)
mmHg	mmHg	Measurement unit the blood pressure(1 mmHg=0.133kPa)
	irregular heartbeat	irregular heartbeat detection
	Current time	Month/Day/Year, Hour/Minute
	Grade	The grade of the blood pressure
	Heartbeat	Heartbeat detection during the measurement

3 things to know about blood pressure

What are systolic pressure and diastolic pressure?

When ventricles contract and pump blood out of the heart, blood pressure reaches its maximum value, the highest pressure in the cycle is known as systolic pressure. When the heart relaxes between heartbeats, the lowest blood pressure is diastolic pressure.


What is the standard blood pressure classification?

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

CAUTION:

Only a physician can tell your normal blood pressure range. Please contact a physician if your measuring result falls out of the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

Blood pressure (mmHg)	Optimal	Normal	High normal	Mild hypertension (Grade I)	Moderate hypertension (Grade II)	Severe hypertension (Grade III)
SYS	< 120	120-129	130-139	140-159	160-179	≥ 180
DIA	< 80	80-84	85-89	90-99	100-109	≥ 110

Irregular heartbeat detector

An irregular heartbeat is detected when a heartbeat rhythm varies while the device is measuring systolic pressure and diastolic pressure. During each measurement, blood pressure monitor will keep a record of all the pulse intervals and calculate the average value of them. If there are two or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 25\%$, or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of $\pm 15\%$, then the irregular heartbeat symbol will appear on the display with the measurement result.

CAUTION:

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

4 before first use

Before you use the appliance for the first time, please do as follows: carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bags and cardboard) out of reach of children. Check the appliance after unpacking for any damage, possible from transportation.


- Open the battery cover.
- Install the batteries by matching the correct polarity, as shown.
Always use the correct battery type (4 x AAA batteries).
- Replace the cover.

Replace the batteries whenever the below happens:

The  +Lo shows

The display dims

The display does not light up

CAUTION:

- Remove batteries if the device is not likely to be used for some time.
- The old batteries are harmful to the environment; do not dispose with other daily trash.
- Remove the old batteries from the device and follow your local recycling guidelines.
- Do not dispose of batteries in fire. Batteries may explode or leak.
- Do not use new and used batteries together.
- Do not use different types of batteries together.

Measurement principle

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a “zero pressure” equivalent to the air pressure. Then it starts inflating the arm cuff and, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsation, which is used to determine the systolic and diastolic pressure, and also the pulse rate.

Date, time and unit settings

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory. (year: 2014 –2054, time format: 24 H/12 H)


1. When the monitor is off, hold pressing “SET” button for about 3 seconds to enter the mode for setting the year. The blinking numeral represents [YEAR].


2. Press the “MEM” button to change the [YEAR].


3. When you get the right year, press “SET” to confirm your selection and it will turn to the next step.


4. Repeat steps 2 and 3 to set [MONTH] and [DAY].

5. Repeat steps 2 and 3 to set [HOUR] and [MINUTE].

6. Repeat steps 2 and 3 to set the [UNIT]. The standard settings are mmHg.



7. After the unit is set, the LCD will display "done" first and then display all the settings you have done before turning it off.

5 measuring blood pressure


Tie the cuff

1. Make your arm clear of clothes.
2. Put the cuff on your upper arm, then position the tube off-centre towards the inner side of the arm in line with the little finger.
3. The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.
4. Position the cuff with the edge 2 to 3 cm above your elbow as indicated on the cuff.
5. Sit comfortably with your tested arm resting on a flat surface.
6. Patients with Hypertension: The middle of the cuff should be at the level of the right atrium of the heart;
Before starting measurement, please sit comfortably with legs uncrossed, feet flat on the floor, back and arm supported.
 - Rest for 5 minutes before measuring.
 - Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
 - The patient must relax as much as possible and do not talk during the measurement procedure.
 - For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

Start the measurement

1. When the monitor is off, press the “START/STOP” button to turn on the monitor, and it will finish the whole measurement automatically.

2. Press the “START/STOP” button to power off; otherwise it will turn off automatically within 1 minute.

6 evaluating results

Irregular heartbeat indication

This instrument can identify possible cardiac arrhythmia disorders (irregular heartbeat) during measurement and if necessary indicates the measurement with the icon . This may be an indicator of arrhythmia. Arrhythmia is a condition where the heart rhythm is abnormal as a result of defects in the bioelectrical system controlling the heartbeat. The symptoms (omitted or premature heart beats, slow or excessively fast heart rate) may be caused, among other things, by heart disease, age, physical predisposition, excessive use of stimulants, stress or lack of sleep. Arrhythmia can only be ascertained through examination by your doctor. Repeat the measurement if the flashing icon is displayed after the measurement. Please note that you should rest for 3 minutes between measurements and not talk or move during the measurement. If the icon appears often, please contact your doctor. Any self-diagnosis and treatment based on the results may be dangerous. It is vital to follow your doctor's instructions.

WHO classification

According to the WHO Guidelines/Definitions and the latest findings, the test results can be classified and evaluated according to the following chart. Your physician will tell you your normal blood pressure range as well as the point at which you will be considered at risk.


Table for classification of blood pressure values (unit of measurement mmHg) for adults:

Blood pressure (mmHg)	Optimal	Normal	High normal	Mild hypertension (Grade I)	Moderate hypertension (Grade II)	Severe hypertension (Grade III)
SYS	< 120	120-129	130-139	140-159	160-179	≥ 180
DIA	< 80	80-84	85-89	90-99	100-109	≥ 110
Action	-	-	Check with doctor	Check with doctor	Seek medical advice	Seek medical advice



7 saved records

After finishing a measurement the appliance will save the result automatically. The max. number of saved measurements is 60. If this is exceeded, it will automatically delete the oldest measurement so that the new measurement can be saved. To recall the saved records, follow the steps below.


- When the monitor is off, press the “MEM” button to show the average value of the latest three measurement records.

- Press the “MEM” button to go up through the saved measurements, or press the “SET” button to go down to the saved measurements.

The record number, date and time will display alternately.

8 delete all saved records

You can delete all records by following the steps below.

1. Hold down the "MEM" button for 3 seconds when the monitor is in the memory recall mode, the flash display will show "dEL ALL".

2. Press "SET" to confirm deleting, the LCD displays "dEL dONE" and the monitor will turn off.


3. If you don't want to delete the records, press "START/STOP" to switch off the appliance without deleting.

4. If all records are deleted successfully, the display will show "___"

Please note that with this function all records will be deleted. It is not possible to delete only one, or a few records.

9 tips for measurement

Measurements may be inaccurate if taken in the following circumstances.

- Results can be inaccurate when measuring within 1 hour after dinner or drinking.
- Results can be inaccurate when measuring immediate after drinking tea or coffee or after smoking.
- Results can be inaccurate when measuring within 20 minutes after taking a bath.
- Results can be inaccurate when measuring when talking or moving your fingers.
- Results can be inaccurate when measuring in a very cold environment.
- Results can be inaccurate when measuring when you want to discharge urine.

10 maintenance

- Clean your blood pressure monitor carefully only with a slightly moistened cloth.
- Do not use detergents or solvents.
- Never immerse the appliance in water or any other liquid to protect against electrical hazards and damage.
- When storing the device, make sure that no heavy objects are placed on top of it. Remove the batteries. The cuff should not have any sharp kinks.
- Always store the blood pressure monitor in a dry place and avoid sunshine.
- Avoid intense shaking and collisions.
- Avoid dusty and unstable temperature environments.
- The cuff is not washable and cannot be cleaned with any liquid.

11 frequently asked questions

Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times every day. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.
2. If the person takes medicine, the pressure will vary more.
3. Wait at least 3 minutes for another measurement.
4. Blood pressure might fluctuate after eating and drinking.

Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise etc., Also, there is the “white coat” effect, which means blood pressure usually increases in clinical settings.

Is the result the same if measuring on the right upper arm?

It is all right for both arms, but there will be some different results for different people. We suggest you measure the same arm every time.

12 troubleshooting

This section includes a list of error messages and frequently asked questions for problems you may encounter with your wrist blood pressure monitor. If the product is not operating as you think it should, check here before arranging for servicing.

Problem	Symptom	Check this	Remedy
No power	Display doesn't show anything.	Batteries are exhausted.	Replace with new batteries.
		Batteries are inserted incorrectly.	Insert the batteries correctly according to the symbols.
Low batteries	 +Lo shown on the display	Batteries are low.	Replace with new batteries
Error message	E1 shows	The inflation is too low.	Refasten the cuff and then measure again.
	E2 shows	The cuff is not secure.	Refasten the cuff and then measure again.
	E3 shows	Can't inflate the air during the process of measurement	Refasten the cuff and then measure again.
	E10 or E11 shows	The monitor detected motion while measuring.	Movement can affect the measurement. Relax for a moment and then measure again.
Error message	E20 shows	The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then measure again.
	E21 shows	Measure incorrectly.	Relax for a moment and then measure again.
	EExx shows on the display	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service for further assistance. Refer to the warranty for contact information and return instructions.
Warning message	"out" shows	Out of measurement range	Relax for a moment. Refasten the cuff and then measure again. If the problem persist, contact your physician.

13 specifications

Power supply	Battery powered mode: 4 x AAA batteries (6V DC)
Display mode	Digital LCD V.A. 40.5 mm x 60 mm
Measurement mode	Oscillographic testing mode
Measurement principle	Oscillometric bloodpressure and heartbeat measuring with automatically 0-setting for atmospheric pressure at every measurement.
Measurement range	Rated cuff pressure: 0 mmHg-299 mmHg (0kPa ~39.9kPa) Measurement pressure: SYS: 60 mmHg~230 mmHg (8.0 kPa~30.7 kPa) DIA: 40 mmHg~130 mmHg (5.3 kPa~17.3 kPa) Pulse value: (40-199) beat/minute
Accuracy	Pressure: 5°C-40°C within ± 0.4 kPa(3 mmHg) Pulse value: $\pm 5\%$
Normal working condition	A temperature range of: +5 °C to +40 °C A relative humidity range of 15% to 90% non-condensing, but not requiring a water vapour partial pressure greater than 50 hPa. An atmospheric pressure range of: 700 hPa tot 1060 hPa.
Storage & transportation condition	Temperature: -20°C to +60°C A relative humidity range of $\leq 93\%$ RH, non-condensing, at a water vapour pressure up to 50hPa
Measurement perimeter of the arm	About 22 cm-32 cm
Net weight	Approx. 168 g (excluding the dry cells)

External dimensions	Approx. 110 mm x 110 mm x 41 mm (exclude the cuff)
Attachment	4 x AAA batteries, user manual and storage box
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP21 (It means the device could be protected against solid foreign objects of 12,5mm Φ and greater, and against vertically falling water drops.)
Software version	A01
Device classification	Internally Powered ME Equipment

14 complied European standards list

Risk management	EN ISO 14971: 2012 Medical devices - Application of risk management to medical devices (ISO14971: 2007, Corrected version 2007-10-01)
Labelling	EN ISO 15223-1: 2016 Symbols for use in the labelling of medical devices
User manual	EN 1041: 2008 Information supplied by the manufacturer of medical devices
General requirements for safety	<p>EN 60601-1: 2006+A1: 2013 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance</p> <p>EN 60601-1-11: 2015/IEC 6060-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment</p> <p>IEC 80601-2-30:2009 + A1:2013 Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers</p>

Electromagnetic compatibility	EN 60601-1-2:2015/ IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests
Performance requirements	EN ISO 81060-1: 2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type EN 1060-3: 1997+A2: 2009 Non-invasive blood pressure Part 3: Supplementary requirements for electromechanical blood pressure measuring system
Clinical investigation	EN 1060-4: 2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers
Usability	EN 60601-1-6:2010+A1:2015/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment -Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability IEC 62366-1:2015 Medical devices - Application of usability engineering to medical devices
Software life-cycle processes	EN 62304: 2006/AC: 2008 / IEC 62304: 2006+ A1: 2015 Medical device software - Software life cycle processes

15 EMC guidance

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.”

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment BDA432, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
2. Guidance and manufacturer's declaration –electromagnetic emissions and Immunity.

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity		
Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact, ±2kV, ±4 kV, ±8kV, ±15 kV air	±8 kV contact, ±2kV, ±4 kV, ±8kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	not applicable	not applicable
Surge IEC61000-4-5	not applicable	not applicable
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	not applicable	not applicable
Magnetic field IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz
Conducted RF IEC61000-4-6	not applicable	not applicable
Radiated RF IEC61000-4-3	10 V/m 80 MHz - 2,7 GHz 80% AM op 1kHz	10 V/m 80 MHz - 2,7 GHz 80% AM op 1kHz
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 3

Guidance and manufacturer's declaration - electromagnetic immunity							
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	385	380-390	TETRA 400	Pulse modulation b) 18Hz	1.8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM c) ± 5kHz deviation 1kHz sine	2	0.3	28
	710	704-787	LTE Band 13,17	Pulse modulation b) 217Hz	0.2	0.3	9
	745						
	780						
	810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation b) 18Hz	2	0.3	28
	870						
	930						
	1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS	Pulse modulation b) 217Hz	2	0.3	28
	1845						
	1970						

Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Frequency (MHz)	Band a) (MHz)	Service a)	Modulation b)	Modulation b) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	2	0.3	28
	5240	5100-5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	0.2	0.3	9
	5500						
	5785						

Indications for use

The blood pressure monitor is digital monitors intended for use in measuring blood pressure and heartbeat rate with arm circumference ranging from 22 cm to 32 cm (about 8¾~12½").

It is intended for adult indoor use only.

Contraindications

- 1.The device should not be used by any person who may be suspected of,or is pregnant .
- 2.The device is not suitable for use on patients with implanted,electrical devices, such as cardiac pacemakers, defibrillators.

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general terms and conditions of service and warranty

We do not need to remind you of the importance of service. After all, we develop our products to a standard so that you can enjoy them for many years, without any concerns. If, nevertheless, there is a problem, we believe you are entitled to a solution straight away. Hence our products come with an exchange service, on top of the rights and claims you are entitled to by law. By exchanging a product or part, we save you time, effort and costs.

2-year full manufacturer's warranty

1. Customers enjoy a 2-year full manufacturer's warranty on all Inventum products. Within this period, a faulty product or part will always be exchanged for a new model, free of charge. In order to claim under the 2-year full manufacturer's warranty, you can either return the product to the shop you bought it from or contact the Inventum customer service department via the form at www.inventum.eu/service-aanvraag.
2. The 2-year warranty period starts from the date the product is bought.
3. In order to claim under the warranty, you must produce a copy of the original receipt.
4. The warranty applies only to normal domestic use of Inventum products within the Netherlands.

Breakdowns or faults outside the warranty period

1. Breakdowns or faults in small or large domestic appliances outside the warranty period, can be reported to the customer services department via the form at www.inventum.eu/service-aanvraag or by calling the customer services department.
2. The customer services department may ask you to send the product for inspection or repair. The costs of dispatch will be at your expense.
3. The inspection to establish whether repair is possible involves a charge. You need to grant your permission for this, in advance.
4. In the event of a large domestic appliance, Inventum, at your request, can send out a service engineer. In that case, you will be charged the call-out costs, as well as parts and labour.

5. In the event of an instruction to repair, the repair costs must be paid in advance. In the event of a repair by a service engineer, the costs of the repair must be settled with the engineer onsite, preferably by means of PIN payment.

Warranty exclusions

1. The following is excluded from the aforesaid warranties:
 - normal wear and tear;
 - improper use or misuse;
 - insufficient maintenance;
 - failure to comply with the operating and maintenance instructions;
 - unprofessional installation or repairs by third parties or the customer himself;
 - non-original parts used by the customer;
 - use for commercial or business purposes;
 - removal of the serial number and/or rating label.
2. In addition, the warranty does not apply to normal consumer goods, such as:
 - dough hooks, baking tins, (carbon) filters, etc.;
 - batteries, bulbs, carbon filters, fat filters etc.;
 - external connection cables;
 - glass accessories and glass parts such as oven doors;
 - and similar items.
3. Transport damage not caused by Inventum is also excluded. Therefore, inspect your new device before starting to use it. If you detect any damage, you must report this to the store where you purchased the product within 5 working days, or to the Inventum customer service department via the form at www.inventum.eu/service-aanvraag. If transport damage is not reported within this period, Inventum does not accept any liability in this respect.
4. The following are excluded from warranty and/or replacement: faults, loss of and damage to the device as a result of an event that is normally insured under the home contents insurance.

Important to know

1. The replacement or repair of a faulty product, or a part thereof, does not lead to an extension of the original warranty period.
2. If a complaint is unfounded, all costs arising from it will be at the customer's expense.
3. Following expiry of the warranty period, all costs of repair or replacement, including administration costs, dispatch and call-out charges, will be charged to the customer.
4. Inventum cannot be held liable for damage caused by external events, unless this liability arises from mandatory statutory provisions.
5. These warranty and service provisions are governed by Dutch law. Disputes will be settled exclusively by the competent Dutch court.

Guangdong Transtek Medical Electronics Co., Ltd.

Zone A, No. 105, Dongli Road, Torch Development District,
528437 Zhongshan, Guangdong, China.

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CE 0123 Model no. TMB-1491-S

Medical Device Safety Service GmbH
Schiffgraben 41, 30175 Hannover, Germany

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Apparaten B.V.**