PERSONENWEGSCHAAL - PERSONENWAAGE - PÈSE PERSONNE - SCALE

INVENTUM

SINCE 1908

PW 601 BM

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ENGLISH

Legislation requires that all electrical and electronic equipment must be collected for reuse and recycling. Electrical and electronic equipment marked with the symbol indicating separate collection of such equipment must be returned to a municipal waste collection point.

The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired.

DESCRIPTION OF THE APPLIANCE

PAGE 20

1. SAFETY INSTRUCTIONS



PAGE 20

2. INTERESTING FACTS

PAGE 20

3. BEFORE FIRST USE
4. BATTERIES
5. WEIGHT MEASUREMENT
6. SETTING USER DATA
7. ACTIVITY LEVELS
8. TAKING MEASUREMENTS
9. DISPLAYING MEASURED VALUE
10. EVALUATION OF RESULTS
11. TIPS
12. INCORRECT MEASUREMENTS
13. CONVERT FROM KG/ST/LB
14. CLEANING & MAINTENANCE

- PAGE 21
PAGE 21
PAGE 21
PAGE 21
PAGE 21
PAGE 22
PAGE 22
PAGE 22
PAGE 23
PAGE 24
PAGE 24
PAGE 24
PAGE 26

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

1. RVS platform
2. Display
3. ON-button
4. Down-key
5. Up-key
6. SET-key (confirm)
7. Battery compartment (underneath)

Minimum/ Maximum capacity: 10 kg to 150 kg

/ ¼ St to 24 St/ 20 lb to 330 lbs.

Division per 100 gram

LCD display: size 40 mm

2 interchangeable batteries, type CR2032 - 3 Volt

1 | SAFETY INSTRUCTIONS

- Please read the instructions carefully and completely before operating this appliance and retain for future use.
- The appliance is intended solely for domestic use within the home.
- Always use batteries which corresponds with the voltage printed on the bottom of the appliance.
- **Attention:** Do not use during pregnancy. Inaccurate measurements may be caused by the amniotic fluid.
- **Attention:** The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired. 
- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Stand on the scale without moving and with equal weight distribution on both feet!
- Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur.
- **Caution!** Do not stand on the scale with wet feet and do not step on the scale when its surface is wet - danger of slipping!
- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).
- When you want to store the scale for a longer period, we advise to keep it in the original carton.
- Always ensure that the scale is placed on a flat, stable surface and only use it indoors.
- Repairs may only be carried out by Inventum customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- If you decide never to use the appliance again, we advise you to dispose of it in special containers or at a

municipal waste collection point.

- **Never immerse the scale in water or other liquids.**
- The scale can accept a maximum load of 150 kg (330 lb, 24St). The weight measurement is displayed in increments of 100 g (0.2lb, ¼St). The measuring results for the body fat, body water and muscle percentages are displayed in 0.1% increments.
- On delivery the scale is set to the units "cm" and "kg". A switch located on the back of the scale enables you to change to "inch", "pound" and "stones" (lb, St).
- The device is only for private use and is not intended for medical or commercial purposes.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- If the appliance is used for other purposes than intended, or is it not handled in accordance with the instruction manual, the full responsibility for any consequences will rest with the user. Any damages to the product or other things are not covered by the warranty.

2 | INTERESTING FACTS

The measuring principle of the diagnostic scale

This scale determines the body weight, bodyfat, body water and the percentage of muscle according to the principle of Bioelectric Impedance Analysis (BIA). With this process, an imperceptible, completely safe and harmless current enables the determination of body tissue proportions within seconds. This measurement of electrical resistance (impedance) and the inclusion of constants and/or individual values (age, height, sex, activity levels) can be used to determine the percentage of fat, water and muscle in the body. Muscle tissue and water have good electrical conductivity, as the fat cells hardly conduct the current due to their very high resistance.

Weighing correctly

If possible, always weigh yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results. Only the long-term trend counts. Stand on the scale without moving and with equal weight distribution on both feet!

Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur. Do not stand on the scale with wet feet and do not step on the scale when its surface is wet - danger of slipping! Place the scale on a firm, flat surface; a firm floor covering is a basic requirement for correct weight measurement.

3 | BEFORE FIRST USE

Before you use the appliance for the first time, please do as follows: carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bags and cardboard) out of reach of children. Check the appliance after unpacking for any damage, possible from transportation. If present, pull the battery insulating strip off the battery compartment cover or remove the battery's protective film and insert the battery according to the polarity.


4 | BATTERIES

If provided, pull out the battery insulating strip at the battery compartment or remove the protective film from the battery and insert the battery with the correct polarity. If the scale fails to operate, remove the battery completely and reinsert it. Your scale is provided with a "Change battery" indicator. If the scale is operating with a battery which is too weak, "Lo" appears on the display and the scale will switch off automatically. In this case, the batteries must be exchanged (2 x 3 V, CR2032). The used, completely drained standard and rechargeable batteries must be disposed of in specially marked collection containers, at toxic waste collection points or electrical product retailers. You are legally obligated to dispose of the batteries.

Note: Batteries containing pollutants are marked with the following symbols: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

5 | WEIGHT MEASUREMENT

Switch the scale on by pressing the "ON"-button. As a self-test, the complete display will appear up to "0.0". The scale is now ready to measure your weight.

Step on the scale and stand still with your weight distributed evenly on both legs. The scale will start measuring your weight immediately and a short time later will display your weight.


When you step off the scale, the scale will switch off after a few seconds.

6 | SETTING USER DATA

To measure your body fat percentage and other physical data, you must enter your personal data.

The scale has 10 memory positions in which you and other members of your family can save and recall personal settings.

Switch on the scale. Wait until the display shows "0.0". Then press "SET". The first memory position will appear flashing on the display.

Now you can enter the following setting:

Memory position - Person	1 to 10
Body size	100 to 220 cm (3'-03" to 7'-03")
Age	10 to 100 years
Sex	male (♂), female (♀)
Degree of activity - Activity	1 to 5

You can enter the relevant settings with short or long pressure on the "up" ▲ button or "down" ▼ button. Confirm the settings in each case by pressing "SET". Then the scale is ready for measuring. If you do not carry out a measurement, the scale switches off automatically after several seconds.

7 | ACTIVITY LEVELS

When selecting the level of activity the average and long-term counts.

- **A 1:** No physical activity.
- **A 2:** Little physical activity. Little and light physical exercises (e.g. walking, light gardening and gymnastics). This mode should also be selected before the start of a training or dietary program. After approx. 6-10 weeks you can switch to the next level of activity.
- **A 3:** Moderate physical activity. 30 minutes of physical exercise 2 to 4 times a week (e.g. jogging, cycling, tennis, etc.). This mode should be selected as soon as you notice a general improvement in your well-being/fitness in the course of a training/dietary program. After 8-12 weeks of continuous training/diet, you can change to the next level of activity.
- **A 4:** Substantial physical activity. 30 minutes of physical exercise 4 to 6 times a week.
- **A 5:** Intensive physical activity. Daily intensive training or physical work (e.g. intensive running, construction work, etc.). For this mode you should have an overall high level of mobility, endurance and power and maintain this level over a longer period of time.

8 | TAKING MEASUREMENTS

After all data have been entered, you can now measure your weight, body fat and the other data.


- Switch on the scale.
- Select the memory preset where your basic personal data are stored by pressing the “up” ▲ button repeatedly. These data are then shown consecutively with the display of the ideal weight until the display “0.0” appears.

- Step on the scale with bare feet and stand still on the electrodes.
- **Important:** There may be no contact between feet, legs, calves and thighs. Otherwise the measurement will not be performed correctly.

The following data are automatically displayed consecutively:

1. Weight and percentage of body fat
2. Weight and percentage of body water


3. Weight and percentage of muscle
4. Bone mass and active metabolic rate “Kcal”

Then your weight and the percentage of body fat are shown again until the display “0.0” appears. The scale is now ready for another measurement (repeat measurement). If you do not carry out a measurement, the scale switches off automatically after several seconds.

Your measured value is automatically stored during switch-off. Two memories are available for each user for this purpose.

9 | DISPLAYING MEASURED VALUE

After selecting a memory and carrying out a determination of body fat, the display “0.0” appears briefly. Pressing the “down” ▼ button displays your last stored measured value. The weight, percentage of body fat, percentage of body water and muscle mass are displayed consecutively. Pressing the “down” ▼ button again displays your second -to-the-last stored measured value. The “M” stands for Memory.

10 | EVALUATION OF RESULTS

Body fat percentage

The following body fat percentages are for guidance (contact your physician for further information). This diagnostic scale is equipped with an interpretation pyramid for simple classification of your body fat values. This shows you at a glance which category your percentage of body fat is classified in.

Man	top	good	average	poor
Age	very good	good	average	poor
10-14	< 11%	11 - 16%	16,1 - 23%	> 23,1%
15-19	< 12%	12 - 17%	17,1 - 22%	> 22,1%
20-29	< 13%	13 - 18%	18,1 - 23%	> 23,1%
30-39	< 14%	14 - 19%	19,1 - 24%	> 24,1%
40-49	< 15%	15 - 20%	20,1 - 25%	> 25,1%
50-59	< 16%	16 - 21%	21,1 - 26%	> 26,1%
60-69	< 17%	17 - 22%	22,1 - 27%	> 27,1%
70-100	< 18%	18 - 23%	23,1 - 28%	> 28,1%

Female	top	good	average	poor
Age	very good	good	average	poor
10-14	< 16%	16 - 21%	21,1 - 26%	> 26,1%
15-19	< 17%	17 - 22%	22,1 - 27%	> 27,1%
20-29	< 18%	18 - 23%	23,1 - 28%	> 28,1%
30-39	< 19%	19 - 24%	24,1 - 29%	> 29,1%
40-49	< 20%	20 - 25%	25,1 - 30%	> 30,1%
50-59	< 21%	21 - 26%	26,1 - 31%	> 31,1%
60-69	< 22%	22 - 27%	27,1 - 32%	> 32,1%
70-100	< 23%	23 - 28%	28,1 - 33%	> 33,1%

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values states. It should, however, be noted that there could be a danger to health in the case of extremely low values.

Body water

The body water percentage is normally within the following ranges:

Man

Age	poor	good	very good
10-100	< 50%	50-65%	> 65%

Woman

Age	poor	good	very good
10-100	< 45%	45-60%	> 60%

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage.

Body water measurement with this scale is not suitable for drawing medical conclusions, for example concerning age-related water retention. If necessary ask your physician. Basically, a high body water percentage should be the aim.

Muscle percentage

The muscle percentage is normally within the following ranges:

Man

Age	low	normal	high
10-14	< 44%	44 - 57%	> 57%
15-19	< 43%	43 - 56%	> 56%
20-29	< 42%	42 - 54%	> 54%
30-39	< 41%	41 - 52%	> 52%
40-49	< 40%	40 - 50%	> 50%
50-59	< 39%	39 - 48%	> 48%
60-69	< 38%	38 - 47%	> 47%
70-100	< 37%	37 - 46%	> 46%

Female

Age	low	normal	high
10-14	< 36%	36 - 43%	> 43%
15-19	< 35%	35 - 41%	> 41%
20-29	< 34%	34 - 39%	> 39%
30-39	< 33%	33 - 38%	> 38%
40-49	< 31%	31 - 36%	> 36%
50-59	< 29%	29 - 34%	> 34%
60-69	< 28%	28 - 33%	> 33%
70-100	< 27%	27 - 32%	> 32%

Bone mass

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise. With appropriate muscle building, you can also strengthen your bone structure.

Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water).

Attention: Please do not confuse bone mass with bone density. Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone

hardness (e.g. osteoporosis) using this scale. Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender).

AMR = Active Metabolic Rate

The energy consumption of a human being rises with increasing physical activity and is measured on the diagnostic scale in relation to the degree of activity entered (1-5). To maintain your existing weight, the amount of energy used must be reintroduced into the body in the form of food and drink. If less energy is introduced than is used over a longer period of time, your body will obtain the difference largely from the amount of fat stored and your weight will reduce. If, on the other hand, over a longer period of time more energy is introduced than the total active metabolic rate (AMR) calculated, your body will be unable to burn off the excess energy, and the excess will be stored in the body as fat and your weight will increase.

Results in relation to time

Note that it is only the long-term trend which counts.

Short period fluctuations in weight over a few days are mostly the result of a loss of fluid.

The interpretation of the results will depend on changes in your:

overall weight and body fat, body water and muscle percentages, as well as on the period during which these changes take place. Rapid changes within days must be distinguished from medium term changes (over weeks) and long term changes (months).

A basic rule is that short term changes in weight almost exclusively represent changes in water content, whereas medium and long term changes may also involve the fat and muscle percentages.

- If your weight reduces over the short term, but your body fat percentage increases or remains the same, you have merely lost water - e.g. after a training session, sauna session or a diet restricted only to rapid weight loss.
- If your weight increases over the medium term and the body fat percentage falls or stays the same, then you could have built up valuable muscle mass.
- If your weight and body fat percentage fall simultaneously then your diet is working - you are losing fat mass. Ideally you should support your diet with physical activity, fitness or power training. By this means, you can increase your muscle percentage over the medium term.

11 | TIPS

Important when measuring body fat/body water/muscle percentage:

- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Stand still during the measurement.
- Wait several hours (6-8) after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.

The measurement is not reliable for:

- Children under approx. 10 years of age;
- Professional sporters and bodybuilders;
- Pregnant women;
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis;
- Persons taking cardiovascular medication. Persons taking vasodilating or vasoconstricting medications;
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

12 | INCORRECT MEASUREMENTS

If the scale detects an error during weighing, “FFFF” or “Err” appears in the display. If you step onto the scale before “0.0” appears in the display, the scale will not operate properly.

Possible causes of errors	Remedy
The maximum load-bearing capacity of 150 kg (330 lbs) was exceeded.	Only weigh the maximum permissible weight.
The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
Your body fat lies outside the measurable range (less than 3% or greater than 55%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.
Your percentage of water lies outside the measurable range (less than 25% or greater than 75%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.

And remember: Physical exercise is the basis for a healthy body.

13 | CONVERT FROM KG/ST/LB

This scale can display measurements in kilograms, Stone and pounds. You can set this with the switch on the bottom of the scale. Kg (kilograms), St (Stone) or Lb (pounds).

14 | CLEANING & MAINTENANCE

Wipe the scale only with a damp cloth. **Never immerse the scale in water or any other liquids.** Handle the scale with care, protect the scale from impact with hard object and do not drop it.

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GUARANTEE & AFTER SALES SERVICE

- This product is guaranteed for a period of 2 years from the date of purchase against mechanical and electrical defect.
- Guarantee can only be given if proof of purchase (e.g. till receipt) can be shown. Therefore keep it in a safe place. The guarantee period starts at date of purchase.
- Guarantee is not valid in case of:
 - normal wear;
 - damage through misuse use or when the instructions were not followed correctly;
 - lack of maintenance;
 - not following the operation- and safety instructions;
 - incompetent montage or repairs made by yourself or others;
 - using accessories and spare parts which are not suitable for the appliance; for assessment by our technical department.
- The guarantee is not valid for:
 - normal wear;
 - external cables;
 - (interchangeable) signal- and controllamps for assessment by our technical department.

Return the appliance, including proof of purchase and description of the complaint, stamped to:

Inventum Holland BV
Generatorstraat 17
3903 LH VEENENDAAL
The Netherlands

or to: new E.T.P. bvba
Passtraat 121
9100 Sint Niklaas
BELGIUM

- Pack it carefully (preferably in the original carton). Ensure unit is clean.
- Enclose your name and address and quote model number on all correspondence.

When the guarantee is no longer valid or when the period of guarantee has passed, we will charge any costs that flow out of an assessment and/or repairing.

On request we will supply a statement before repairing. Without requesting a statement, we assume that you, by sending the appliance, automatically order to repair the appliance.

A guarantee period of two months is valid for repairs done by our technical department and the used spare parts, counting from the repair date. When a malfunction occurs within this period which directly flows out of the repairs and the used spare parts, the consumer has the right of repairs without any costs.

Inventum Holland BV
Generatorstraat 17, 3903 LH , Veenendaal, Nederland
Tel.: 0318-542222, Fax.: 0318-542022, Mail: helpdesk@inventum.eu

Inventum Belgium BVBA
Xavier de Cocklaan 66 bus 8, 9830, Deurle, België
Tel: 09 220 30 20, Fax.: 09 220 30 22, Mail: info@inventumbelgium.be

Website: www.inventum.eu