INVENTUM

maakt 't moment

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English

The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired.

1. safety instructions page 11
 2. description of the appliance page 13
 3. interesting facts page 13
 4. before first use page 14
 5. personal data page 14
 6. weigh only page 15
 7. body analysis page 15
 8. evaluation of results page 15
 9. tips page 17
 10. incorrect measurements page 17
 11. Cleaning & maintenance page 17
- Guarantee & after sales Service* page 33

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1 safety instructions


- **Please read the instructions carefully and completely before operating this appliance and retain for future use.**
- The appliance is intended solely for domestic use within the home.
- Always use batteries which corresponds with the voltage printed on the bottom of the appliance.
- **Attention:** Do not use during pregnancy. Inaccurate measurements may be caused by the amniotic fluid.
-  **Attention:** The scale may not be used by persons with medical implants (e.g. pacemakers). Otherwise their function may be impaired.
- Caution! Do not stand on the outer edge to one side of the scale otherwise it may tip!
- An electrical appliance is not a toy and should therefore always be placed outside of the reach of children. Children are not always able to understand potential risks. Teach children a responsible handling of electrical appliances.
- Children should be supervised to ensure that they do not play with the appliance.
- This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- If possible, always weigh yourself at the same time of day (preferably in the morning), after using the toilet, on an empty stomach and unclothed, in order to achieve comparable results.
- Stand on the scale without moving and with equal weight distribution on both feet!
- Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur.
- Protect the scale against hard knocks, moisture, dust, chemicals, major temperature fluctuations and heat sources which are too close (stove, heating radiators).
- Always ensure that the scale is placed on a flat, stable surface and only use it indoors.
- If you decide never to use the appliance again, we advise you to dispose of it in special containers or at a municipal waste collection

point.

- Never immerse the scale in water or other liquids.
- The scale can accept a maximum load of 180 kg (397 lb, 28St: 5lb). The weight measurement is displayed in increments of 100 g (0.2lb, 1 lb). The measuring results for the body fat, body water and muscle percentages are displayed in 0.1% increments.
- On delivery the scale is set to the units “cm” and “kg”. For changing the measuring unit, please read the information in chapter BEFORE FIRST USE.
- When you want to store the scale for a longer period, we advise to keep it in the original carton.
- Repairs may only be carried out by Inventum customer service or authorized dealers. Before submitting any complaint, first check the batteries and replace them if necessary.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- The device is only for private use and is not intended for medical or commercial purposes.
- Please note that technically related measuring tolerances are possible, as the scale is not calibrated for professional medical applications.
- If the appliance is used for other purposes than intended, or is it not handled in accordance with the instruction manual, the full responsibility for any consequences will rest with the user. Any damages to the product or other things are not covered by the warranty.
-  Legislation requires that all electrical and electronic equipment must be collected for reuse and recycling. Electrical and electronic equipment marked with the symbol indicating separate collection of such equipment must be returned to a municipal waste collection point.

2 appliance description

1. Display
2. Keys for setting **SET** ▲ ▼
3. Glass surface
4. Measuring electrodes
5. Battery compartment (bottom)

5

3 interesting facts

The measuring principle of the diagnostic scale

This scale determines the body weight, bodyfat, body water and the percentage of muscle according to the principle of Bioelectric Impedance Analysis (BIA). With this process, an imperceptible, completely safe and harmless current enables the determination of body tissue proportions within seconds. This measurement of electrical resistance (impedance) and the inclusion of constants and/or individual values (age, height, gender) can be used to determine the percentage of fat, water and muscle in the body. Muscle tissue and water have good electrical conductivity, as the fat cells hardly conduct the current due to their very high resistance.

Weighing correctly

If possible, always weigh yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results. Minor deviations are normal. When you drink 1 glass of water for example a difference of 300 grams can occur. Place the scale on a firm, flat surface; a firm floor covering is a basic requirement for correct weight measurement.

Your scale is an accurate instrument. The human body is the hardest thing to weigh, because it is always in motion. Therefore a 99% accuracy of the body weight is the highest possible. Stand on the scale without moving and with equal weight distribution on both feet!

Do not underestimate the weight of clothing and when you have just eaten or drunken anything. Only the long-term trend counts.

4 before first use

Before you use the appliance for the first time, please do as follows: carefully unpack the appliance and remove all the packaging material. Keep the material (plastic bags and cardboard) out of reach of children. Check the appliance after unpacking for any damage, possible from transportation. Place the scale on a firm, flat surface. Remove the batteries from the protective packaging and insert with the correct polarity. See the diagram inside the battery compartment.

Press the UNIT button on the back of the scale to select the correct measuring unit (kg/st/lb).

kg = kilograms
lb = pounds
st = stone

Place the scale on a flat and stable surface and wait for the scale to turn itself off.

Before the first use the scale must be calibrated. This must be done every time you use the scale on a different spot or when the batterie is replaced.

Place the scale on a flat and stable surface. With your foot press on the middle of the glass surface. The full display will show and when it shows "0.0" the scale is ready for use.

5 personal data


You must enter personal user parameters in order to measure your body fat percentage and other physical data. The scale has 10 memory positions in which you and other members of your family can save and retrieve personal settings. Now you can enter the following setting:

Memory position (P)	P0 to P9
Body size	100 to 240 cm/3 ft 3 inch - 7 ft 10inch
Age	10 to 80 years
Gender	male (♂), female (♀)
Sportsman or -woman	sportman (♂) or sportswoman (♀) For people who exercise regularly.


- Press the SET key when the display is switched off to change the setting.

- Press the SET key again to select the correct user, press the ▲ or ▼ keys and press the SET key again to confirm.

- After you have entered the correct user, the symbol for gender will blink. To select the gender (male, female, sportsman, sportswoman), press the ▲ or ▼ keys and press the SET key again to confirm.

- After you have entered the gender, the indication for body size (height) will blink (preset 165 cm). To set your own height press the ▲ or ▼ keys and press the SET key again to confirm.

- After you have entered your body size, the indication for age will blink (preset 30). To set the correct age (male, female, sportsman, sportswoman), press the ▲ or ▼ keys and press the SET key again to confirm.

After all parameters have been entered, you can now measure your weight, body fat and the other data. The display will show 0,0 kg.

Stand on the scale barefoot immediately and let the scale perform a measuring. When you step on the scale the next time the scale will recognize you. When you step off the scale, the scale will switch off after a few seconds.

TIP: To reset a memory position for another user perform the same steps as described above. All positions are preset.

6 weigh only

The scale has a step-on function. Step barefoot onto the scale. Stand still with your weight distributed evenly on both legs. The scale will start measuring your weight immediately and a short time later will display your weight. When you step off the scale it will turn off after a few seconds.

7 body analysis


When you stay barefoot on the scale after it displayed your weight, the scale will automatically perform a body analysis. The scale can identify users according to the weight.

Important: When using the analysis function, make sure there is no contact between your feet, legs, calves or thighs. Otherwise the measurement will not be performed correctly. After the weight is stable, the "O" signal moves from left to right for two times, and then shows the data.

The scale is measuring the bodyfat, -water, etc. The scale can identify a person by their weight. If it detects more users with a similar weight, it will display the relative user positions. Use the ▲ or ▼ keys to choose the right user. Stay on the scale to complete the measurement.

The scale will now display the results: Bodyfat, bodywater, BMI, muscle, bone mass, calories and weight.


To give you the opportunity to see all the results, the scale will show them two times before it turns off.

It will also display the user information again for confirmation.

The scale can identify a person by their weight. If it detects more users with a similar weight (< 2kg), it will display the relative user positions. Use the ▲ or ▼ keys to choose the right user.

Stay on the scale to complete the measurement.

8 evaluation of the results

The percentage of bodyfat represents the in terms of percentage part of the complete body weight. There is no real advise for this value, because it depends on gender and age. Bodyfat is very important for your body, to ensure good functioning. It protects organs, hinges, regulates the body temperature, stores vitamins and acts as a reserve for the body. Females need a bodyfat percentage between 10% and 13%. Males need a bodyfat percentage between 2% and 5%, with exception of sports people.

The following body fat percentages are guide values (contact your physician for further information).

Male

age	low	normal	high	very high
10-39	< 8	8 - 20	20 - 25	> 25
40-59	< 11	11 - 22	22 - 28	> 28
60-80	< 13	13 - 25	25 - 30	> 30

Female

age	low	normal	high	very high
10-39	< 21	21 - 33	33 - 39	> 39
40-59	< 23	23 - 34	34 - 40	> 40
60-80	< 24	24 - 36	36 - 42	> 42

A lower value is often found in athletes. Depending on the type of sports, training intensity and physical constitution, values may result which are below the recommended values states. It should, however, be noted that there could be a danger to health in the case of extremely low values.

A large part of the human body contains water. This bodywater can be found through the entire body. About 60% of the bodymass consists of water with men. With women this is about 55% (caused by a higher bodyfat percentage).

The body water percentage is normally within the following ranges:

Male

age	poor	good	very good
10 - 15	< 58	58 - 72	> 72
16 - 30	< 53	53 - 67	> 67
31 - 60	< 47	47 - 61	> 61
61 - 80	< 42	42 - 56	> 56

Female

age	poor	good	very good
10 - 15	< 57	57 - 67	> 67
16 - 30	< 47	47 - 57	> 57
31 - 60	< 42	42 - 52	> 52
61 - 80	< 37	37 - 47	> 47

Body fat contains relatively little water. Therefore persons with a high body fat percentage have body water percentages below the recommended values. With endurance athletes, however, the recommended values could be exceeded due to low fat percentages and high muscle percentage. This scale is unsuitable for measuring body water in order to draw medical conclusions concerning age-related water retention, for example. If necessary ask your physician.

The BMI (Body Mass Index) is normally within the following ranges:

under weight	normal	over weight	obese
< 18.5	18.5 - 25	25 - 30	> 30

The muscle percentage is normally within the following ranges:

age	male healthy	female healthy
10-80	> 40	> 34

Like the rest of our body, our bones are subject to the natural development, degeneration and ageing processes. Bone mass increases rapidly in childhood and reaches its maximum between 30 and 40 years of age. Bone mass reduces slightly with increasing age. You can reduce this degeneration somewhat with healthy nutrition (particularly calcium and vitamin D) and regular exercise. With appropriate muscle

building, you can also strengthen your bone structure. Note that this scale will not show you the calcium content of your bones, but will measure the weight of all bone constituents (organic substances, inorganic substances and water).

The bone mass is normally within the following ranges:

	male			female		
weight	< 60 kg	60-75 kg	> 75 kg	< 45 kg	45-60 kg	> 60 kg
bonemass	2,5 kg	2,9 kg	3,2 kg	1,8 kg	2,2 kg	2,5 kg

Attention: Please do not confuse bone mass with bone density. Bone density can be determined only by means of a medical examination (e.g. computer tomography, ultrasound). It is therefore not possible to draw conclusions concerning changes to the bones and bone hardness (e.g. osteoporosis) using this scale.

Little influence can be exerted on bone mass, but it will vary slightly within the influencing factors (weight, height, age, gender).

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The basal metabolic rate is the amount of energy required by the body at complete rest to maintain its basic functions. This value largely depends on weight, height and age. This is displayed on the scale as the amount of energy that your body requires in any case and it must be reintroduced into your body in the form of nutrition. If you take on less energy over the longer term, this can be harmful to your health.

CALORIES ANALYSIS (Unit: kcal)					
	Age	< 17	18 - 29	30 - 50	> 50
MALE	< 50 kg	1350	1200	1115	1075
	50 - 65 kg	1485	1320	1227	1183
	65 - 80 kg	1890	1680	1561	1505
	80 - 95 kg	2295	2040	1896	1828
	> 95 kg	2700	2400	2230	2150
FEMALE	< 50 kg	1265	1180	1085	1035
	50 - 65 kg	1392	1298	1194	1139
	65 - 80 kg	1771	1652	1519	1449
	80 - 95 kg	2151	2006	1845	1760
	> 95 kg	2530	2360	2170	2070

9 tips

Remember that only long-term trends are important. Short-term fluctuations in weight over a few days are usually the result of a loss of fluid.

Important when measuring body fat/body water/muscle percentage:

- If possible, always weighs yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results.
- The measurement may only be carried out while barefoot and it is helpful if the soles of your feet are slightly damp. Completely dry soles can result in unsatisfactory results, as they have insufficient conductivity.
- Stand still during the measurement.
- Wait several hours after unusually strenuous activity.
- Wait approx. 15 minutes after getting out of bed so that the water in your body can be distributed.
- Remember that only long-term trends are important. Short-term changes in weight within a few days are usually caused by a loss of fluid. Body water plays an important role towards our general wellbeing.

The measurement is not reliable for:

- Children under approx. 10 years of age;
- Professional sporters and bodybuilders;
- Pregnant women;
- Persons with fever, undergoing dialysis, with symptoms of edema or osteoporosis;
- Persons taking cardiovascular medication.
- Persons taking vasodilating or vasoconstricting medications;
- Persons with substantial anatomical deviations in the legs relative to their total height (leg length considerably shorter or longer than usual).

10 incorrect measurement

If the scale detects an error during measurement, this is indicated with **Err**.

Possible causes of errors	Causes of error
The electrical resistance between the electrodes and the soles of your feet is too high (e.g. with heavily callused skin).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary. Remove the calluses on the soles of your feet if necessary.
The fat percentage is outside the measurable range (less than 5% or greater than 50%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.
The water percentage is outside the measurable range (less than 10% or greater than 80%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.

Possible causes of errors	Causes of error
The muscle and bone percentages are outside the measurable range - age and gender related (less than 10% or greater than 80%).	Repeat weighing barefoot. Slightly moisten the soles of your feet if necessary.
The maximum load-bearing capacity of 180 kg (397 lb) was exceeded [Err].	Only weigh the maximum permissible weight.
You are not standing correctly on the scale [Err2].	Make sure your feet are placed on the measuring strips.

11 cleaning & maintenance

The scale should be cleaned occasionally. Dirt and dust has a tendency to collect on the rubber electrodes. Always stand on the scale barefoot or wearing socks as the rubber surfaces will quickly become dirty if you wear shoes. Use a damp cloth to clean the scale. Do not use abrasive detergents and never immerse the scale in water.

- Do not use aggressive solvents or cleaning agents!
- Never immerse the scale in water or any other liquids.
- Handle the scale with care, protect the scale from impact with hard object and do not drop it.

Your scale is fitted with a change battery indicator. If the scale is operating with batteries which are too weak, Lo appears on the display and the scale will switch off automatically. In this case, the batteries must be exchanged (2 x CR2032, 3V). If necessary please use a small screwdriver to remove the batteries. Batteries can contain toxins that are harmful to health and the environment. The used, completely drained standard and rechargeable batteries must be disposed of in specially marked collection containers, at toxic waste collection points or electrical product retailers. You are legally obligated to dispose of the batteries.

Note: Batteries containing pollutants are marked with the following symbols: Pb = Battery contains lead, Cd = Battery contains cadmium, Hg = Battery contains mercury.

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general terms and conditions of service and warranty

We do not need to remind you of the importance of service. After all, we develop our products to a standard so that you can enjoy them for many years, without any concerns. If, nevertheless, there is a problem, we believe you are entitled to a solution straight away. Hence our products come with an exchange service, on top of the rights and claims you are entitled to by law. By exchanging a product or part, we save you time, effort and costs.

2-year full manufacturer's warranty

1. Customers enjoy a 2-year full manufacturer's warranty on all Inventum products. Within this period, a faulty product or part will always be exchanged for a new model, free of charge. In order to claim under the 2-year full manufacturer's warranty, you can either return the product to the shop you bought it from or contact the Inventum customer service department via the form at www.inventum.eu/service-aanvraag.
2. The 2-year warranty period starts from the date the product is bought.
3. In order to claim under the warranty, you must produce a copy of the original receipt.
4. The warranty applies only to normal domestic use of Inventum products within the Netherlands.

Breakdowns or faults outside the warranty period

1. Breakdowns or faults in small or large domestic appliances outside the warranty period, can be reported to the customer services department via the form at www.inventum.eu/service-aanvraag or by calling the customer services department.
2. The customer services department may ask you to send the product for inspection or repair. The costs of dispatch will be at your expense.
3. The inspection to establish whether repair is possible involves a charge. You need to grant your permission for this, in advance.
4. In the event of a large domestic appliance, Inventum, at your request, can send out a service engineer. In that case, you will be charged the call-out costs, as well as parts and labour.
5. In the event of an instruction to repair, the repair costs must be paid in advance. In the event of a repair by a service engineer, the costs of the repair must be settled with the engineer onsite, preferably by means of PIN payment.

Warranty exclusions

1. The following is excluded from the aforesaid warranties:
 - normal wear and tear;
 - improper use or misuse;
 - insufficient maintenance;
 - failure to comply with the operating and maintenance instructions;
 - unprofessional installation or repairs by third parties or the customer himself;
 - non-original parts used by the customer;
 - use for commercial or business purposes;
 - removal of the serial number and/or rating label.
2. In addition, the warranty does not apply to normal consumer goods, such as:
 - dough hooks, baking tins, (carbon) filters, etc.;
 - batteries, bulbs, carbon filters, fat filters etc.;
 - external connection cables;
 - glass accessories and glass parts such as oven doors;
 - and similar items.
3. Transport damage not caused by Inventum is also excluded. Therefore, inspect your new device before starting to use it. If you detect any damage, you must report this to the store where you purchased the product within 5 working days, or to the Inventum customer service department via the form at www.inventum.eu/service-aanvraag. If transport damage is not reported within this period, Inventum does not accept any liability in this respect.
4. The following are excluded from warranty and/or replacement: faults, loss of and damage to the device as a result of an event that is normally insured under the home contents insurance.

Important to know

1. The replacement or repair of a faulty product, or a part thereof, does not lead to an extension of the original warranty period.
2. If a complaint is unfounded, all costs arising from it will be at the customer's expense.
3. Following expiry of the warranty period, all costs of repair or replacement, including administration costs, dispatch and call-out charges, will be charged to the customer.
4. Inventum cannot be held liable for damage caused by external events, unless this liability arises from mandatory statutory provisions.
5. These warranty and service provisions are governed by Dutch law. Disputes will be settled exclusively by the competent Dutch court.

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Apparaten B.V.**